





WoundFix™

Amniotic Membrane Skin Substitute
for Outpatient Wound Care

Q4217




Treat Your Patients with the Advantage of Natural Tissue Engineering®

WoundFix™ is a FastActing® dehydrated amniotic membrane skin substitute.

Complete viable BioECM® preserved by the HydraTek® process provides all the native placental components the human body needs for its natural ability to:

✓ Support Correct Tissue Formation

WoundFix™ DFU / VLU Study Highlights



98%
closure achieved
at 12 weeks

1.5
sq. cm

Average
Diabetic
Ulcer Size

7.4
sq. cm

Average
Wound Size
Closed with
WoundFix™

Average Number of
WoundFix™ Applications

2.5

The HydraTek® Advantage

The Preserved Source.

Maintains natural properties providing a market-leading allograft transplant.

Collagen Content Retained

Unprocessed Placental Tissue




441.6 µ/mg

HydraTek® Placental Tissue



435 µ/mg

Deviation



WoundFix™

- Amnion-only membrane
- Thin & versatile
- Ultra FastActing®
- Omni-directional application
- Room temp storage

The HydraTek® Advantage



Allograft Potential

Retains high levels of native properties & collagens



Immune Privileged



FastActing®

Designed to resorb faster & go to work quicker



Antimicrobial

Placental tissues have been shown to be bacteriostatic against certain bacteria in vitro

Ordering Information

| WoundFix™ | Size | WoundFix™ | Size |
|-----------|----------|-----------|-----------|
| WF001 | 1x1cm | WF096 | 96 sq cm |
| WF004 | 4 sq cm | WF128 | 128 sq cm |
| WF008 | 8 sq cm | WF160 | 160 sq cm |
| WF016 | 16 sq cm | WF192 | 192 sq cm |
| WF024 | 24 sq cm | WF243 | 243 sq cm |
| WF032 | 32 sq cm | WF324 | 324 sq cm |
| WF064 | 64 sq cm | WF405 | 405 sq cm |
| WF081 | 81 sq cm | WF486 | 486 sq cm |